

## HOW TO USE

The OPTI® CCA-TS2 analyzer is easy to run:

- The **user-friendly touch screen**, with simple menu navigation, can be used with minimal training.
- **Automated sample aspiration** reduces staff time and training.
- **Heated measuring chamber** allows for testing of iced samples.

## Running the OPTI CCA-TS2 analyzer:



### 1. Scan the cassette bar code.



2. Remove the cassette from its packaging.



3. Place the cassette into the sample measurement chamber and close the cover.



4. Aspirate the sample. When sample aspiration and measurement is complete, the results will automatically appear on the screen.